

## **DRÄGER Babylog 2, short-term ventilation device for premature babies and newborn infants ...**

Article number: 1001535

**780,00 €**  
net price



### **Product description**

DRÄGER Babylog 2, short-term ventilation device for premature babies and newborn infants to carry out controlled ventilation (IPPV) with PEEP or spontaneous respiration with CPAP during transportation or for resuscitation, inspiration time-controlled or triggered manually, pressure-limited expiration is time-controlled, inspiration pressure 10-60 mbar plateau pressure, expiration pressure 0-10 mbar, ventilatory rate 10-60/min., driving gas O<sub>2</sub> or dried, oil-free compressed air 2-6 mbar, H:195 mm x W:120 mm x D:200 mm, mit DRÄGER Babymix 2 O<sub>2</sub> mixer, second-hand

-> [Click here to come to the site.](#) <-

