

DRÄGER Oxylog 2000 transport ventilator with accessories, time-controlled, constant ...

Article number: 1001559

2.450,00 €
net price



Product description

DRÄGER Oxylog 2000 transport ventilator with accessories, time-controlled, constant volume for patients with a tidal volume from 1000 ml ventilation mode: IPPV \ SIMV \ CPAP. Variables 1:E, also IRV, with integrated PEEP valve 0-15 mbar. Air-Mix \ No-Air-Mix switch, integrated expiratory volume measurement with display, second-hand

-> [Click here to come to the site.](#) <-

