

**MELAG Euroklav 23VS+, chamber size  
45 cm x diameter 25 cm, 22 liter, load  
capacity 4 kg ...**

Article number: 1003751

**3.650,00 €**  
net price



**Product description**

MELAG Euroklav 23VS+, chamber size 45 cm x diameter 25 cm, 22 liter, load capacity 4 kg of instruments, 1 kg textiles, dimensions (WxHxD): 42.5 x 48.5 x 69 cm, the unit fits on a 60 cm tabletop, weight 45 kg, power consumption 2600 Watt, used

-> [Click here to come to the site.](#) <-

