

## Wisap CO2 Aqua-Purator 1611. The device provides continuous flushing of the abdomen and ...

Article number: 1002131

**480,00 €**  
net price



### Product description

Wisap CO2 Aqua-Purator 1611. The device provides continuous flushing of the abdomen and aqua dissection during a diagnostic laparoscopy. The accessory (bottle pin for NACL) is not included in the delivery and must be ordered directly from the manufacturer. Technical data: continuously adjustable pressure up to 750 mmHg. NACL instillation quantity: selectable 1 litre or Hi-Flo 3 l/min, second-hand

-> [Click here to come to the site.](#) <-

