

**OLYMPUS OTV-SI: video processor  
with integrated cold light source 150  
Watt halogen ...**

Article number: 1003939

**2.900,00 €**  
net price



**Product description**

OLYMPUS OTV-SI: video processor with integrated cold light source 150 Watt halogen (MD-151), for connection of video endoscopes: ENF-V2, ENF-VT2, ENF-VQ, or of camera heads: OTV-S7H-N, OTV-S7H-1N. OTV-S7H-1D, used

-> [Click here to come to the site.](#) <-

