

Olympus PSD-10 VET. Monopolar HF device, the handle can be activated using the foot ...

Article number: 1003148

780,00 €
net price



Product description

Olympus PSD-10 VET. Monopolar HF device, the handle can be activated using the foot switch, CUT 80 W, blend 55 W, coag 40 W, used. THIS HF SURGICAL DEVICE IS SUPPLIED WITH ACCESSORIES (NEUTRAL PLATE CABLE FOR SINGLE NEUTRAL ELECTRODES, MONOPOLAR ELECTRODE, 3 DIFFERENT INSERTS and FOOTSWITCH) This item is exclusively designed for being used in the field of veterinary medicine or for export purposes

-> [Click here to come to the site.](#) <-

