

DRÄGER Photo therapy 800, with chassis, for the light irradiation of premature or newly ...

Article number: 1001357

590,00 €
net price



Product description

DRÄGER Photo therapy 800, with chassis, for the light irradiation of premature or newly born infants in the incubator or in the infrared-heated thermal bed to reduce the serum-bilirubin concentration (phototherapy), frequency: 50 Hz, output: 230W, current: 2.2 A, lamp housing 54 x 25.5 x 20.5cm, second-hand

-> [Click here to come to the site.](#) <-

