

DRÄGER 5400 transport incubator on transport cradle, for external and internal O2 and ...

Article number: 1002974

3.200,00 €
net price



Product description

DRÄGER 5400 transport incubator on transport cradle, for external and internal O2 and air supply, bottle connections for O2 + air, power supply via mains connection or car battery, with Babylog 2000, dimensions without cradle height: 86 cm x width: 150cm x depth: 56cm, incl. bronchial suction (O2 injector), dimensions of cradle H 86 x W 189 x D 61cm, very good condition, used

-> [Click here to come to the site.](#) <-

