

BRAUN Perfusor Compact, very good condition, delivery rate 0.1-99.9 ml/h, bolus rate 800 ...

Article number: 1001335

690,00 €
net price



Product description

BRAUN Perfusor Compact, very good condition, delivery rate 0.1-99.9 ml/h, bolus rate 800 ml/h, new battery, second-hand
This perfusor is ready for use with B. BRAUN perfusor syringes. If you want to use different syringe types let us know please!
For retroactive changement (after sale) we have to charge EUR 55,-

-> [Click here to come to the site.](#) <-

