

BRAUN Perfusor fm, very good condition, delivery rate 0.1 to 200 ml/hour, with new ...

Article number: 1001337

490,00 €
net price



Product description

BRAUN Perfusor fm, very good condition, delivery rate 0.1 to 200 ml/hour, with new battery and integrated mains component, second-hand This perfusor is ready for use with B. BRAUN perfusor syringes. If you want to use different syringe types let us know please! For retroactive chngement (after sale) we have to charge EUR 55,-

-> [Click here to come to the site.](#) <-

