

**BRAUN Perfusor Compact S, volume:
0.01 - 200 ml/hour, time selection
0.01-99h, max. ...**

Article number: 1002673

750,00 €
net price



Product description

BRAUN Perfusor Compact S, volume: 0.01 - 200 ml/hour, time selection 0.01-99h, max. infusion rate: 200 ml/h, max. bolus rate: 1200 ml/h, new battery, used This perfusor is ready for use with B. BRAUN perfusor syringes. If you want to use different syringe types let us know please! For retroactive changement (after sale) we have to charge EUR 55,-

-> [Click here to come to the site.](#) <-

