

**CAREFUSION Alaris plus GH, syringe pump, compatible with many standard disposable ...**

Article number: 1005117

**690,00 €**  
net price



**Product description**

CAREFUSION Alaris plus GH, syringe pump, compatible with many standard disposable syringes 50ml (Braun Omnifix, BD, Codan, Fresenius, Terumo). Max. Infusion rate with 50ml. Syringes 0.1ml/h-1200ml/h. Max. Bolus rate 10ml/h-1200ml/h. System accuracy +/-2%. Rechargeable NIMH battery to operate the device for approx. 4 hours, used

-> [Click here to come to the site.](#) <-

