

**FRESENIUS Injectomat TIVA Agilia.
Syringe pump for battery and mains
operation, incl. ...**

Article number: 1005291

980,00 €
net price



Product description

FRESENIUS Injectomat TIVA Agilia. Syringe pump for battery and mains operation, incl. holding claw for infusion stand. Approx. 50 usable syringe types. Delivery rate 0.1 to 99.9ml/h in 0.1ml increments. Bolus rate 50-1200ml. Please note that the medication database must be programmed by Fresenius field service. second hand

-> [Click here to come to the site.](#) <-

