

**FRESENIUS Injectomat TIVA Agilia.  
Syringe pump for battery and mains  
operation, incl. ...**

Article number: 1005291

**980,00 €**  
net price



**Product description**

FRESENIUS Injectomat TIVA Agilia. Syringe pump for battery and mains operation, incl. holding claw for infusion stand. Approx. 50 usable syringe types. Delivery rate 0.1 to 99.9ml/h in 0.1ml increments. Bolus rate 50-1200ml. Please note that the medication database must be programmed by Fresenius field service. second hand

-> [Click here to come to the site.](#) <-

