

**SMITH EQ-5000 Level 1 Equator, convective patient warmer. Adjustable temperature range: ...**

Article number: 1005201

**590,00 €**  
net price



**Product description**

SMITH EQ-5000 Level 1 Equator, convective patient warmer. Adjustable temperature range: 36°C, 40°C, 44°C. Delivery incl. 2-metre supply hose, second-hand. We recommend the use of convective warming pads from Moeck&Moeck. These are compatible with almost all patient warming systems. Moeck&Moeck also offers special warming pads for veterinary medicine. Please request a quotation directly from Moeck&Moeck: [vet@moeckundmoeck.de](mailto:vet@moeckundmoeck.de)

-> [Click here to come to the site.](#) <-

